



Eye Drop Vaccination

Foreye drop vaccination treatment of day old poultry chicks



Droplet Stay Function (droplet does not fall)



Confirmation by Checking

Credence Biospray Vaccinator



Credence Eye Dropper Vaccinator

Eye Dropper vaccinator, is a very simple and portable method of administering live vaccines directly into the eyes day old poultry chicks. The System also count the number of chicks vaccinated.

Equipment Features

- Long Lasting User Friendly Equipment
- Sturdy metallic steel frame
- Push Button / Sensor operated
- Electronic Adjustable dosage
- Accurate and easy to set dosage from 5 ml to 30 ml
- Electrical system ensures continues operation.
- Counting function of Total. No. of Birds
- Cooling Chamber for Live Vaccine
- Portable Equipment

How it Works

- Whenever the operator presses the foot switch or if the sensor sensed the Device will deliver a drop and keep it suspended. **Note :** The drop will not fall and will keep suspending on the nozzle tube.

- When the operator positions the eye of the chick near the nozzle the droplet is absorbed in the eye of the chick.
- Also the count will go up by 1 value.
- The Vaccine temperature is maintained at cool temp. with the help of gel filled jackets.

Shipping Setup

- Chick Eye Dropper Equipment
- Consumable Spares kit , Tool Kit



Check out the Eye Dropper video @ Youtube

Technical Specifications

Capacity	- 1,800 chicks per hour (approx.)
Target	- Day old chicks
Technique	- Eye Drop Vaccination
Power Requirement	- 230 V, 50 W , Single Phase (optional battery)
Weight	- 12 kg.
Dimensions (LxWxH)	- 500 x 180 mm x 250 mm (body dimensions)
Dosages	- 0.02 ml - 0.05 ml per chick box.
Models	- Electric



Credence Indomechanics India

C102, Shivane , Pune 411023
Maharashtra, India

+917767056868 , +919890293069
credenceindo@gmail.com

www.credencindo.com